



Oxygen Concentrator Safety and Cleaning Guide

Electrical Safety

- Make sure the electrical outlet is working properly
- Make sure the plug is securely plugged in the electrical outlet to avoid overheating and reduce the risk of fire
- Do not use extension cords or adaptors to connect the device
- In case of doubt, have your electrical outlet checked by a professional

Cleaning

Your concentrator consists of a **device (filter and bubble humidifier)** and a **nasal circuit (nasal cannula or mask)**. Unplug your device before cleaning. Here are some tips:

Filter

- Clean the filter every other day. To remove the filter, gently pull on the filter located on the back or the side of the machine
- Wash with mild dish detergent
- Rinse thoroughly and dry with a clean towel
- Wait until completely dry before re-installing

Bubble Humidifier

- Change the water every **other day**. Only use distilled water.
- **Once a week**, remove the bubble humidifier, wash in a solution of warm water and dish detergent, then rinse with warm water
- **Once a month**, soak for 30 minutes in a solution of three parts warm water to one part white vinegar. Rinse thoroughly and fill with distilled water to the fill line as indicated by the manufacturer. Do not overfill.

Nasal Circuit (Nasal Cannula or Mask)

- **Once a day**, wipe mask or nasal cannula with a disinfecting wipe or a damp cloth
- **Once a week**, wash mask or nasal cannula as well as the exterior of the tubing with soapy warm water. Rinse thoroughly and allow to air dry.

OSR Médical Inc.
Biomedical Services

Main Branch

1400 Sauvé Street West, suite 1512
Montréal, QC H4N 1C5
Tel. : 514 332-6387
1 877 433-3956

www.osrmedical.com