



Safety Precautions When Using Oxygen at Home

Oxygen will not explode; however, it will make a fire burn faster. In other words, in the presence of oxygen, a small fire can become dangerous and get out of control.

The following safety precautions must be observed when using an oxygen machine at home:

- Put up a notice saying that you are using oxygen
- Never smoke
- Never use matches or a lighter
- Never use candles or fondue burners
- Never grease or oil the unit
- Never store the oxygen concentrator or tank near a source of heat (heating duct, fire place, wood stove, etc.)
- Never open an oven door or work over a hot stove while using oxygen
- Make sure the tank is securely fastened

NOTE

I understand the importance of these security precautions.
It is my responsibility to share this information with my family and those close to me.

On this day: _____

I hereby acknowledge having read and understood the above safety precautions.

Signed:

(Client)

(OSR Médical Inc. employee)

OSR Médical Inc.
Biomedical Services

Main Branch

1400 Sauvé Street West, suite 1512
Montréal, QC H4N 1C5
Tel. : 514 332-6387
1 877 433-3956

www.osrmedical.com