



Safety Precautions When Using Oxygen at Home

IMPORTANT INFORMATION FOR CLIENTS AND OSR MEDICAL'S TEAM:

Oxygen will not explode; however, it will make a fire burn faster. In other words, in the presence of oxygen, a small fire can become dangerous and get out of control.

The safety precautions below must be followed when using an oxygen machine at home:

- Put a notice saying that you are using oxygen
- **Never** smoke
- Never use matches or a lighter
- Never use candles or fondue burners
- Never grease or oil the unit
- Never put oily or greasy creams on your face and on your hands
- Never store the oxygen concentrator or tank near a source of heat (heating duct, fireplace, wood stove, etc.)
- Never open an oven door or work over a hot stove while using oxygen
- Make sure the tank is securely fastened
- Never use adhesive tape on valves or on the oxygen tank
- Make sure to follow the information provided by OSR Medical's clinician and/or technician before manipulating the oxygen tank and any other related equipment
- Keep the oxygen equipment clean
- Use gloves or make sure to have clean non greasy hands when you assemble the oxygen equipment
- Wear clean non greasy clothes
- Open slowly the valves of the oxygen tank

You know the importance of the back system in case of an emergency

NOTE:

I understand the importance of these security precautions

It is my responsibility to share this information with my family and those close to me

I hereby acknowledge having read and understood the above safety precautions

Name: _____ Witness: _____

Signature: _____ Signature: _____

Date: ____ / ____ / ____
DD MM YYYY

Date: ____ / ____ / ____
DD MM YYYY